



Operational Procedure Manual			
Doc no	PRL-OPM- F001	Revision 0.0	Effective Date: 1 Sept 2020
Written By: Gugu Ditshego		Reviewed by: Gugu Ditshego	Approved by: Jabulani Kubheka
Date: Aug 2020		Date: Aug 2020	Date: 31 Aug 2020
	18069 Ngungunyane Street Kwa Thema 1575  Tel: 011 737 1045/2518 Fax: 086 604 3815 Email: info@proteuslab.co.za		

### Verification reports for ABO & Rh Grouping, Methods and Equipment

Equipment / Method	ABO & Rh Manual
Performed by	Gugu Ditshego
Date of Verification	20 February 2021
Date of Report	20 February 2021
Reviewed by	Gugu Ditshego 
Approved by	Jabulani Kubheka
Date Approved	28 February 2021

#### Introduction

Validation was done using PRL-OPM-001. Samples used were fresh whole blood drawn in EDTA tubes. Samples were analysed within sample stability. A comparability study was done using 10 samples. 100 % correlation was achieved. Raw data available to verify traceability & validity of the report.

Results: acceptance criteria are based on 100% correlation. These was achieved.

Statistical analysis: 10 samples were obtained for ABO & Rh grouping

Limitations of the method: The amount of sample volume can limit ABO & Rh typing. Weak positives result due to reagent/cell strength can cause limitation of method.

Conclusion: A 100% correlation was achieved. Therefore, the ABO Method used in the lab is fit for purpose.

Literature: The ABO and Rh **blood grouping** system is based on agglutination reaction. When red **blood** cells carrying one or both the antigens are exposed to the corresponding antibodies, they interact with each other to form visible agglutination or clumping.

### ABO and Rh Data

<u>Proteus Lab</u>	<u>Results</u>	<u>Comparative method</u>	<u>Results</u>
148927	O POSITIVE	35420824	O POSITIVE
995601	A POSITIVE	35420825	A POSITIVE
148795	AB POSITIVE	35420826	AB POSITIVE
995597	B POSITIVE	35420823	B POSITIVE
104936	A POSITIVE	93174434	A POSITIVE
149538	B POSITIVE	35421568	B POSITIVE
995761	O NEGATIVE	35421567	O NEGATIVE
995720	A POSITIVE	35421566	A POSITIVE
149508	O POSITIVE	35421565	O POSITIVE
149560	A NEGATIVE	35421564	A NEGATIVE
<b>Achieved outcome was 100% Correlation.</b>			

**Statistical processing**

Malaria field staining technique	Diagnostic accuracy criteria		Total
	Known Positive	Known Negative	
<b>Positive</b>	# true positive – TP	#false positive – FP	TP+FP
<b>Negative</b>	#false negative – FN	# true negative - TN	FN+TN
<b>Total</b>	TP+FN	FP+TN	N

$$\text{Sensitivity} = \frac{\text{no of true positives}}{\text{no of true positives + no of false negatives}} \times 100$$

$$\text{Specificity} = \frac{\text{no of true negatives}}{\text{no of true negatives + no of false positives}} \times 100$$

$\text{Sensitivity} = \frac{TP}{(TP+FN)} \times 100$ $= \frac{10}{(10+0)} \times 100$ $= \frac{10}{10} \times 100$ $= 100\%$	$\text{Specificity} = \frac{TN}{(TN+FP)} \times 100$ $= \frac{10}{(10+0)} \times 100$ $= \frac{10}{10} \times 100$ $= 100\%$
--	--

		Gold standard		Total	
		Positive	Negative		
The Test Result	Positive	True Positive (TP) 10	False Positive (FP) 0	10	<b>Sensitivity</b> 100%
	Negative	False Negative (FN) 0	True Negative (TN) 10	10	<b>Specificity</b> 100%
	Total	10	10	20	

