

EQA/IQC:

PROBLEM STATEMENT: CK QC > 2SD

POSSIBLE CAUSES (YES/ NO)

REAGENTS EXPIRED	NO
CALIBRATION ACCEPTABLE	NO
IQC ACCEPTABLE	NO
MACHINE SERVICED/ MAINTENANCE PERFORMED	YES
TECHNICAL STAFF COMPETENT	YES

COMMENT:

Recalibrated the analyte and re-ran QC
and it was acceptable. NO patient were
processed.

WAS CORRECTIVE ACTION APPLIED EFFECTIVE: YES/NO

SIGNED OFF BY: Mr

DATE: 24/02/21

Index 02/24/2021 15:33
Print Time [02/24/2021 16:05]

Proteus Lab QC&CAL

Operator [MAINTENANCE

] Device No. 2121907

Page. 6

S.No. Q001 S.ID
Type Serum
QC No. 7

HDL - 2 1.71 : , , , | LDL - 2 4.37 : , , , |

S.No. Q001 S.ID
Type Serum
QC No. 8

CO2 - 1 12.3 : , , , |

S.No. Q001 S.ID
Type Serum
QC No. 9

CO2 - 2 22.7 : , , , |

S.No. Q001 S.ID
Type Serum
QC No. 10

CO2 - 3 30.1 : , , , |

S.No. Q001 S.ID
Type Urine
QC No. 21

CRE-E - 1 4.942 : , , , | U-ALb - 1 11.60 : , , , |

S.No. Q001 S.ID
Type Urine
QC No. 22

CRE-E - 2 9.772 : , , , | U-ALb - 2 130.30 : , , , |

Index 02/24/2021 15:33
Print Time [02/24/2021 16:05]

Proteus Lab QC&CAL

Operator [MAINTENANCE

] Device No. 2121907 Page. 9

S.No. Q002 S.ID
Type Serum
QC No. 1

CK - 1 160.5 : , , , |

S.No. Q002 S.ID
Type Serum
QC No. 2

CK - 2 375.4 : , , , |

Index 02/24/2021 17:02
Print Time [02/24/2021 18:18]

Proteus Lab QC&CAL

Operator [MAINTENANCE

] Device No. 2121907

Page. 3

S.No. Q002 S.ID
Type Urine
QC No. 21

CRE-E - 1 5.475 : , , | U-Alb - 1 13.68 : , , |

S.No. Q002 S.ID
Type Urine
QC No. 22

CRE-E - 2 10.765 : , , | U-Alb - 2 143.35 : , , |

Index 02/24/2021 17:02
Print Time [02/24/2021 17:21]

Proteus Lab QC&CAL

Operator [MAINTENANCE

] Device No. 2121907 Page. 3

S.No. Q001 S.ID
Type Serum
QC No. 1

BUN - 1 6.46 : , , , |

S.No. Q001 S.ID
Type Serum
QC No. 2

BUN - 2 27.95 : , , , |